


Integration PCR: PCR Strategy to Confirm the Genomic Integration Site of Gene Trap Cassettes

An alternative, more direct way to confirm the correct genomic integration site of any ES cell line in our collection is to use the sequence which we initially retrieved by iPCR and used to map the exact genomic location of the mutagen. By extracting the whole genomic sequence/locus for a particular gene from e.g. MGI, you can easily design primers* surrounding the "mapped sequence" and integration site, respectively. By doing so, you can run 4 PCR reactions on gDNA extracted from wildtype/unmodified and mutated ESCs in the following way: by using the primers flanking the potential genomic integration site on wildtype gDNA (wt, black and purple arrows; purple lane in A), you should obtain a single PCR band, confirming that the primers in fact work. Using the same primer pair on gDNA from the respective clone (mut, black and green arrows; green lane in A), you will not obtain a PCR band as the mutagen is too large for PCR amplification. Conversely, under the assumption that the genomic location of the mutagen was mapped correctly, you should obtain a PCR product with a mutagen-specific primer (mut, black and orange arrows; orange lane in A), but not on wt gDNA (wt, black and red arrows; red lane in A). This is graphically summarized below and a typical example shown (A):



Should you observe any deviation from the expected pattern (e.g. if you receive a PCR product in the green lane, or no PCR product in the orange lane), you should be cautious and compare the results also to the bcPCR.

Possible explanations are (i) the retrieved sequence of the INT site was apparently mapped to the wrong chromosomal location, (ii) the clone is a mixed clone, i.e. it contains another clone, which could have been introduced during the production process, or (iii) the clone was initially not targeted homozygously.

Please make sure that you choose the right combination of primers in dependence of the orientation of the mutagen and gene locus, respectively (see below preferred primer combinations are marked in red).

Mutagen-specific PCR primer sequences for

- Retro, Lenti, Tol2-GT:

RV-I NT-fwd GCCAGAACCAGAAGGAACTTGAC
 RV-I NT-rev TACAGACGCAGGCGCATAACATC
 LV-I NT-rev AGAGCTCCTCTGGTTTCCCTTTC
 T2-I NT-fwd GAGCCAGAACCAGAAGGAACTTG
 • Tol 2-pA
 T2-I NT-rev CCGGGCAATGGATTGATATTGCG
 T2-I NT-fwd GAGCCAGAACCAGAAGGAACTTG

Integration site-specific PCR primer sequences

Int-F, custom ~23 bps

Int-R, custom ~23 bps

designed by Clone Manager 9 Software (Sci-Ed Software)

Primer combinations in dependence of gene trap insertion:



Rel. Orientation Strand of Integration	+	-	+	-
Tol2-pA_3'	Int-F + T2-INT-rev OR: T2-INT-fwd + Int-R	Int-R + T2-INT-rev OR: T2-INT-fwd + Int-F	Int-F + T2-INT-rev OR: T2-INT-fwd + Int-R	Int-R + T2-INT-rev OR: T2-INT-fwd + Int-F
Tol2-GT_3'	T2-INT-fwd + Int-R	T2-INT-fwd + Int-F	T2-INT-fwd + Int-R	T2-INT-fwd + Int-F
Retro_5'	Int-F + RV-INT-fwd OR: RV-INT-rev + Int-R	Int-R + RV-INT-fwd OR: RV-INT-rev + Int-F	Int-F + RV-INT-fwd OR: RV-INT-rev + Int-R	Int-R + RV-INT-fwd OR: RV-INT-rev + Int-F
Lenti_5'	Int-F + RV-INT-fwd OR: LV-INT-rev + Int-R	Int-R + RV-INT-fwd OR: LV-INT-rev + Int-F	Int-F + RV-INT-fwd OR: LV-INT-rev + Int-R	Int-R + RV-INT-fwd OR: LV-INT-rev + Int-F

INT-PCR reaction (4 PCR rxns w/ respective primer combinations on gDNA from wt and mut)

crude ES cell lysate 5 µl
 FWD primer 5 µM 2 µl
 REV primer 5 µM 2 µl
 10 mM dNTP mix 1 µl
 10x Kentaq buffer 5 µl

2

