

**Production of Retrovirus for gene trap flipping:
 Calcium-phosphate transfection of packaging cells and
 production of virus stocks**

Adapted from Johannes Zuber

Materials & Solutions

- retroviral packaging cells (PlatE, Cellbiolabs)
- Cre/Flp-expressing plasmid w/ a FACS-able color and/or a selection marker (e.g. MLP-mCherry-Cre-puro)
- optional: second plasmid as "negative control" (e.g. MLP-GFP-puro)
- RV-Helper plasmid

| | |
|---|--|
| <p>2 x HBS (1000ml) <i>italic substances are optional !</i></p> | <p>280 mM NaCl (16.4g) 50 mM HEPES MW 238.31 (11.92g) 1.5 mM Na₂HPO₄ (0.21g anhydr./ 0.58g x12 H₂O) 12 mM Dextrose (2.0g) or Glucose (2.2g) 10 mM KCl (0.74g) H₂O (ad 1000ml)</p> <p>- distribute into 4 x 250ml and adjust pH to to 6.96, 7.00, 7.04, 7.08 (!!!) with 0.5M NaOH</p> <p>- perform a test precipitation with 15µg GFP+ and look for fine precipitate in hemocytometer</p> <p>- with those yielding a fine precipitate perform independent test transfections</p> <p>- working buffers should be filtered through a 0.2µm filter and stored in aliquots (10 ml) in -20°C</p> <p>- thaw always in room temperature (not in 37°C !!!)</p> |
| <p>1M CaCl₂ (20ml)</p> | <p>. 2.94g CaCl₂ ad 20ml H₂O . filter through a 0.2µm filter . store aliquots at -20°C</p> |

Culture of PlatE packaging cells

- medium: ESCM w/o LiF
- for the first couple passages, culture in ESCM w/o LiF + puro 1:10'000 + blast 1:1'000, switch to medium w/o puro/blast prior to transfection

- culture in 5% CO₂, passage ~1:4 every 2-3 days, cells detach easy from plate!
- always discard cells after >20 passages and avoid long episodes of confluence
- 4h prior to transfection, plate ~19 x 10⁶ in 35mL media per 15cm plate. No further medium change before transfection. Transfect as soon as they are adherent and 75-85% confluent.

Calcium-phosphate transfection of packaging cells

quantities for one 15cm plate:

Solution A (total 1500μl, keep pipetting order):

| | |
|---|-------|
| - plasmid DNA | 60μg |
| - helper DNA | 20μg |
| - H ₂ O (1125 - vol. of DNA) | ...μl |
| - 1M CaCl ₂ | 375μl |

Solution B (total 1500μl): - 2 x HBS 1500μl

- prepare both solutions in separate tubes or flasks
- mix dropwise (A) to (B): Use a vortexer as a means to "blow" bubbles through solution B. Use a pasteur pipet for drop-by-drop adding solution A to solution B while the vortexer "blows" bubbles. The mixture should turn turbid.
- leave mixture at RT for 15min
- if precipitate in the end gets flocculent re-suspend it by pipetting up and down or shearing through a 21G needle
- add precipitate to cells dropwise while moving the plate in circular motions (caution: be gently as the packaging cells easily detach from the plate)
- optional: look for a fine precipitate on the cells 15 min after addition
- 8-14h (i.e. o/n) post transfection change media for recovery and/or 8-12h prior (first) infection change media to target cell media (ESCM for ESCs). Choose the volume needed for infection (e.g. 3ml per target 6well).
- optional: check transfection efficiency under a fluorescence microscope
- virus collection time: 24-72h post transfection; highest titers are at 32-60h post transfection
- carefully aspirate viral supernatant from packaging cells, filter through 0.45μm filter and use for infection of ESCs (pls see protocol 'Gene trap flipping')

Production of virus stock

To always have some virus at hand you can produce a batch and store it at -80deg for approximately half a year.

first day after transfection:

- am: change medium of PlatE to 20mL ESCM for 1st pm harvest
- optional: check transfection efficiency under a fluorescence microscope
- pm (24h post transfection), 1st virus harvest: carefully suck off supernatant and filter through 0.45 μ m filter, store on ice o/n in cold room; add another 20mL ESCM to plates for next harvest

second day after transfection:

- am, 2nd harvest: combine w/ 1st harvest; add ESCM to plates for next harvest
- transfer combined supernatant to centrifuge tubes (appr. 35mL per tube) and concentrate by centrifugation (4h @ 24'000rpm), resuspend viral pellets in 0.5ml ESCM/centrifuge tube and pool, store on ice in cold room
- wash/rinse empty centrifuge tubes w/ PBS and keep in hood o/n for virus inactivation
- pm, 3rd harvest: store on ice o/n in cold room; add ESCM to plates for next harvest

third day after transfection:

- am, 4th harvest: combine w/ 3rd harvest and concentrate as before
- discard PlatE
- aliquot concentrated virus stocks (250 μ L)
- freeze down and store at -80deg

To determine future dilution of the virus, do a test-infection on ES cells with several virus dilutions and analyze infection efficiency by FACS 48h post infection.